

For Mitsui Onnetsuki
MI-II and MI-III

My Hand RAKURAKU

II

Born from idea of an Onnetsuki fan,
My Hand RAKURAKU II can deliver heat
to the target area with ease.

It assists you heat the hard-to-apply area such as back and hips.
The flexibility of RAKURAKU II helps you apply the Onnetsuki to
your body more firmly.



How to attach the Onnetsuki to My Hand RAKURAKU II



1. Hold the electric cord of Onnetsuki like this photo



2. Put the electric cord into the chase



3. Slowly push the head of heating part into RAKURAKU II

Precautions

1. Use Onnetsuki with fabric cover on the head of the heating part
2. Do not put the heating part still on one area for long. Move the area after a few seconds.
3. Move the heating part slowly and slightly.

How to use on the BACK

Thoroughly heat around the spine which is hard to do without RAKURAKU II

- ① Move the heating part to the target area
- ② Push the end of the grip forward slowly (*do not pull it downward)



How to use on the HIPS

Use RAKURAKU II horizontally

- ① Hold RAKURAKU II horizontally and twist it around the torso. Move the heating part to the target area.
- ② Hold RAKURAKU II with one hand and push RAKURAKU II forward with the other hand and twist the torso slightly.



RAKURAKU II can use everyday

My Hand RAKURAKU II is easy to hold.



While working at the desk.



While reading on the chair.



While lying down on the bed.

Safety Precaution

1. This product is designed only for Mitsui Onnetsuki MI-II and MI-III. Please do not use for other purposes.
2. Please do not harm the electric cord of Onnetsuki when you attach My Hand RAKURAKU II.
3. If the electric cord of Onnetsuki is pulled too much, it may harm the Onnetsuki.
4. It may burn your skin or cause injury if you heat one area for too long.
5. Please do not use or store in high temperature place. It may cause damage or fire.
6. Keep out of the reach of children.

Distributed by:
Mitsui Onnetsu Co., Ltd.
3-36-1 Hiyoshidai, Tomisato-city, Chiba, Japan 286-0201
TEL 0476-37-8771